

# TYR TACTICAL® MALE EPIC ASSAULTER'S PLATE CARRIER



## ENHANCED DESIGN FEATURES:

### QUICK RELEASE

- Single-handed 4-wire cutaway

### BREATHABILITY & MOISTURE PROTECTION

- Antimicrobial/FR treated padded spacer mesh system

### SCALABLE SYSTEM

- Compatible with our complete line of armored accessories
- Zippers for Assaulter's Zip-On Back Panels
- MOLLE Assaulter's Zip-On Back Panel
- XFrame® Dynamic Load Carriage® compatible
- PALS webbing for MOLLE attachment
- Semi-rigid, lace adjustable, exterior MOLLE cummerbund

### PATENTED INTEGRATED BALLISTIC FRAMEWORK

- Increases structural integrity & reduces back face deformation

### ELONGATED MATERIAL LIFE & ABRASION RESISTANCE

- Engineered utilizing patented PV® material



### PLATE POCKET STYLE

- Incorporates a Gusseted Plate Pocket which relieves tension and pressure normally placed on the body from hard plates in classic plate pocket configuration
- Available in several different shapes depending on the profile (thickness) of your plate

### ADDITIONAL FEATURES

- EPIC® low profile shoulder pads
- Patented Ballistic Vein®
- Hydration/Communication tabs for routing
- Reinforced drag handle
- Repair buckle kit

NIJ 0101.06 Certified Level IIIA  
TYR-M-EPIC-TXP3A

NIJ 0101.06 Certified Level II  
TYR-M-EPIC-NFA2/10

NIJ 0101.06 Certified Level III  
TYR-M-EPIC-NFA3/10



# TYR TACTICAL® BODY ARMOR SIZING & ARMOR SPEC

**MALE CARRIER SIZE RELATES TO THE PLATE SIZE & OPERATORS CHEST WIDTH:**


<p><b>Small:</b> Chest Width: 33" - 37" Plate Pocket: 8.5" x 11.5" Small Plate</p> <p><b>Medium:</b> Chest Width: 37" - 41" Plate Pocket: 9.3" x 12.3" Medium Plate</p> <p><b>Large:</b> Chest Width: 41" - 45" Plate Pocket: 10" x 13" Large Plate</p>	<p><b>XLarge:</b> Chest Width: 45" - 49" Plate Pocket: 10.9" x 13.8" XLarge Plate</p> <p><b>2XLarge:</b> Chest Width: 49" - 53" Plate Pocket: 10.9" x 13.8" XLarge Plate</p>
---	--



**PLATE POCKET STYLES:**

**Profile 1:** Fit's SAPI/ESAPI style plates  
**Profile 6:** Fit's TYR Tactical® Hard Armor Plates

## TYR TACTICAL® SOFT ARMOR SPECIFICATIONS



SOFT ARMOR PACKAGES	X418	NFA2/10*	NFA3/10*	TXP3A
THICKNESS	.23"	.21"	.24"	.21"
WEIGHT/ AREAL DENSITY	.78 lbs psf	.82 lbs psf	1.09 lbs psf	1.23 lbs psf
TESTING STANDARD	NIJ 0101.06 USMC Testing Standards	NIJ 0101.06 FBI Testing Standards	NIJ 0101.06 FBI Testing Standards	NIJ 0101.06
PROTECTION LEVEL	Level II USMC Testing Requirements	Certified Level II FBI & DEA Requirements USSOCOM Fragmentation Requirements	Certified Level IIIA FBI & DEA Requirements	Certified Level IIIA
V <sup>50</sup>	9mm / 1999 fps .357 Mag / 1852	9mm / 1686 fps .357 Mag / 1741 fps	.357 Sig / 1865 fps .44 Mag / 1690 fps	.357 Sig / 1923 fps .44 Mag / 1647 fps
BACK-FACE DEFORMATION	9mm / 33.7mm .357 Mag / 40.5mm	9mm / 33.8mm .357 Mag / 36.6mm	.357 Sig / 30.3mm .44 Mag / 35.4mm	.357 Sig / 28.3mm .44 Mag / 34.5mm
IDEAL USE	Tactical/ Military	Concealable/ Tactical	Tactical	Tactical

\*Available in Gender Neutral and Female Version

# TYR TACTICAL® BODY ARMOR FEATURES



**ONLY AT TYR TACTICAL® PATENTED INTEGRATED BALLISTIC FRAMEWORK**

- 🌀 **Increases Ballistic Protection**
- 🌀 **Reduces Back Face Deformation**
- 🌀 **Provides Moisture Protection**
- 🌀 **Improves Structural Integrity**
- 🌀 **Increases End-User Ergonomics**
- 🌀 **Enhances Performance**
- 🌀 **Elongates Material Life**
- 🌀 **Exceedingly Improves Abrasion Resistance**



SPACER MESH

PATENTED BALLISTIC VEIN®

BALLISTIC SOFT ARMOR

BALLISTIC PLATE

PATENTED PV®

NAME	FUNCTION
SPACER MESH	The body facing material of our plate carriers utilize advanced antimicrobial spacer mesh which is IR/FR treated.
PATENTED BALLISTIC VEIN®	The patented Ballistic Vein® increases load carriage performance, provides a rigid platform, reduces fatigue, stabilizes and enhances soft armor inserts.
SOFT ARMOR	Helps absorb the impact from firearm projectiles and fragmentations.
HARD PLATE	Provides additional protection from rifle rounds.
PATENTED PV®	Increases ballistic capability of the nylon carrier, provides structure, improves both abrasion and fire resistance.



# TYR TACTICAL® BODY ARMOR CLEANING INSTRUCTION



## BALLISTIC PANELS

Removable ballistic panels should be cleaned in the following method:

- Do not immerse in liquid
- Gently wipe-down ballistic panels with a damp sponge or wash cloth using mild detergent. Remove any excess detergent with a damp sponge or wash cloth using clean water.
- Allow to thoroughly air dry before inserting panels into carrier.
- Do not dry clean, machine wash or bleach
- Do not dry outdoors in the sun




## OUTER CARRIER

The nylon carrier is treated with chemicals to reduce your IR signature. Therefore you should never use a phosphorus based detergent to clean your vest, doing so will breakdown its IR treatment causing it to glow under IR light.

- Remove ballistic panels from the front and back of the carrier.
- Brush off as much surface dust as possible.
- Scrub set-in stains with a nylon brush using a mild non-phosphorus detergent such as Dawn dish soap.
- Air-dry indoors at room temperature.

Additional Instructions:

- Do Not Bleach
- Carrier Must Be Completely Dry Before Inserting The Ballistic Panels



# TYR TACTICAL® EPIC® ARMOR KIT



## SOFT ARMOR OPTIONS:

**NFA2/10**

NIJ 0101.06 Certified Level II

**NFA3/10**

NIJ 0101.06 Certified Level IIIA


**TXP3A**

NIJ 0101.06 Certified Level IIIA

## MAKE IT A KIT



FIGURE	DESCRIPTION	ITEM#
A	EPIC® CUTAWAY THROAT	TYR-F-ETHR-CW
B	EPIC® CUTAWAY COLLAR/YOKE	TYR-F-CLRYK-CW
C	THOR SHOULDERS	TYR-LW107-THOR
D	BALLISTIC BICEPS/DELTOID	TYR-MST115
E	EPIC® ASSAULTER'S PLATE CARRIER	TYR-F-EPIC
F	BALLISTIC LOWER ABDOMEN	TYR-LW109MV
G	MV BALLISTIC GROIN	TYR-LW110MV
H	XFRAME® GUNFIGHTER® BELT	TYR-XF-BKBLB
I	INNER BASE BELT	TYR-DBB008
J	XFRAME® STAYS	TYR-XF-STAYS



\*INCLUDED WITH PLATE CARRIER

**AVAILABLE IN MALE & FEMALE VERSIONS**